Case Report Series

Guided Implant Surgery with Simultaneous
Connective Tissue Grafting in the Premolar Region





Dr. Stijn Kregting, Centrum Implantologie Castricum, Netherlands www.cicastricum.nl



The Baseline Situation

Patient with two missing premolars

A 45-year-old healthy female patient with good dental hygiene presented at the dental office for implant therapy. **Premolars 24 and 25 had been extracted five years earlier** due to endodontic failure and a vertical tooth fracture.






Clinical and radiological assessments confirmed the presence of adequate vertical and horizontal bone volume

The Treatment Strategy

1. Guided surgery

Fully guided insertion of two SPI®ELEMENT RC INICELL® Implants (PF 4.0, L 11 mm) followed by immediate restoration








2. Subepithelial tissue graft

For mucosa thickening and papilla restoration, a connective tissue graft (6x14 mm) harvested from the palate (Regio 14-15) was placed and secured with mattress sutures during initial surgery









3. Monolithic zirconia crowns

Loading of the implants with the final restoration after 4 months of healing





The Outcome

After four months, both implants had successfully osseointegrated with healthy soft tissue around the implants. The final zirconia restorations yielded satisfactory esthetic results





