

# Thommen Medical Library Installation DWOS

---



SWISS CRAFTSMANSHIP FOR DENTAL PROFESSIONALS

**Thommen Medical AG** | Headquarters

Neckarsulmstrasse 28 | CH-2540 Grenchen | Tel. +41 (0)32 644 30 20 | Fax +41 (0)32 644 30 25 |  
info@thommenmedical.com | www.thommenmedical.com


## 1 Instruction for use

Electronic instructions for use for Thommen Medical CAD/CAM Products can be found here:


Thommen Ti-Base:	<a href="https://www.online-ifu.com/THM61123">https://www.online-ifu.com/THM61123</a>
Thommen VARIOflex/temp:	<a href="https://www.online-ifu.com/THM61119">https://www.online-ifu.com/THM61119</a>
Thommen PF 3.0 VARIOunite:	<a href="https://www.online-ifu.com/THM61146">https://www.online-ifu.com/THM61146</a>
Thommen Milling Abutment for CAD/CAM R02:	<a href="https://www.online-ifu.com/THM61137">https://www.online-ifu.com/THM61137</a>
Digital Impression (Scanabutment and Analog):	<a href="https://www.online-ifu.com/THM61143">https://www.online-ifu.com/THM61143</a>
DyTi-Base	<a href="https://www.das-eifu.com/">https://www.das-eifu.com/</a>

## 2 Importing the implant library


- Open the Implant Kit Editor




- Import the new library by clicking «+» and «Import Implant Manufacturer Kits»



- Select the .zip file that contains the library



- After a successful import, the new library can be found in the library list (1). Check if the new library is located in the «Thommen Medical» folder and move it if necessary
- Assign the materials you want to use by right clicking on the abutment kit (2)
- Select «Edit supported materials» (3)
- Checkmark the materials you want to assign (4)
- Repeat for all PF sizes in the abutment kit




## 2 Settings for printed models and the Thommen Medical Analog for CAD/CAM


To achieve optimal results with implant analogs sockets in printed models, DWOS allows for adjusting offset parameters. Offset values may be different, depending on the model printer used. Adjust the offset value to these settings:

Print service provider	Printer	Value
Dreve	Prodways	-0.025
3dmedical Print	Stratasys	0.000
-	Formlabs2	0.000

The setting is found in «User Preferences» (1), «Preferences» (2), «Model Builder Options» (3), «Analog Holder Radial Offset» (4)



In the release DWOS 7.2 the offset will be adjustable in the «Virtual Model Builder» tab of the model material settings.



## 3 Contact

For questions or feedback please contact:  
[patrick.waelchli@thommenmedical.com](mailto:patrick.waelchli@thommenmedical.com)  
+41 32 644 30 29

