

Straight to the goal.




Guided Surgery

PF 3.5 – PF 4.5
L 6.5 – L 12.5

Simply brilliant

Fully guided surgery with Thommen Medical Implants enables you to achieve the planned restoration efficiently, in just a few steps. Benefit from improved treatment comfort for patients and dental health practitioners and simplify the treatment process for the practice team.



The All-rounder

Precision


One system for **CONTACT** and **ELEMENT**

Easy handling

Individual design

One system for surgery and preparation

Space for optional instruments



Start your planning with your preferred planning system:

Company	Software	Fully Guided	Pilot
3Shape	Implant Studio	✓	✓
Swissmeda	SMOP	✓	✓
3DIEMME	RealGuide	✓	✓
SICAT GmbH	GALILEOS Implant & SICAT Implant	✓	✓
Planmeca	Romexis	✓	
360dps™	360imaging	✓	✓
Exocad	Exoplan	✓	✓
Zirkonzahn	Zirkonzahn ImplantPlaner	✓	✓

...

The guided surgery system will be available for other planning systems.

For an up-to-date summary visit www.thommenmedical.com, or contact your country's sales representative.

For further information on use, please refer to operating instructions www.ifu-tm.com/THM61144 «Guided Surgery».

Headquarters

Thommen Medical AG
Neckarsulmstrasse 28
2540 Grenchen | Switzerland
Tel. +41 61 965 90 20
Fax +41 61 965 90 21
info@thommenmedical.com