

ÜBERSICHT BOHRPROTOKOLL CONTACT

PF Ø 3.5 enossal Ø 2.7	PF Ø 4.0 enossal Ø 3.5	PF Ø 4.5 enossal Ø 3.5	PF Ø 5.0 enossal Ø 4.2	PF Ø 6.0 enossal Ø 5.0	Drehzahl U/min
2.0 Pilotbohrer	2.0 Pilotbohrer	2.0 Pilotbohrer	2.0 Pilotbohrer	2.0 Pilotbohrer	800
↓	2.8 VECTOdrill	2.8 VECTOdrill	2.8 VECTOdrill	2.8 VECTOdrill	600
	↓	↓	3.5 VECTOdrill	3.5 VECTOdrill	500
			4.3 VECTOdrill	400	
2.7/3.5 CONTACT Profilbohrer	3.5/4.0 CONTACT Profilbohrer	3.5/4.5 CONTACT Profilbohrer	4.2/5.0 CONTACT Profilbohrer	5.0/6.0 CONTACT Profilbohrer	250-300

