

PROTOCOLE DE FORAGE **CONTACT**

PF Ø 3.5 Ø endo-osseux 2.7	PF Ø 4.0 Ø endo-osseux 3.5	PF Ø 4.5 Ø endo-osseux 3.5	PF Ø 5.0 Ø endo-osseux 4.2	PF Ø 6.0 Ø endo-osseux 5.0	Vitesse de rotation t/min
2.0 Foret pilote	2.0 Foret pilote	2.0 Foret pilote	2.0 Foret pilote	2.0 Foret pilote	800
↓	2.8 VECTOdrill	2.8 VECTOdrill	2.8 VECTOdrill	2.8 VECTOdrill	600
	↓	↓	↓	3.5 VECTOdrill	500
4.3 VECTOdrill				400	
2.7/3.5 Foret de profil CONTACT	3.5/4.0 Foret de profil CONTACT	3.5/4.5 Foret de profil CONTACT	4.2/5.0 Foret de profil CONTACT	5.0/6.0 Foret de profil CONTACT	250-300

